Mo Yag Laser (2 Cases)	<b>Office</b>	/ Date	<b>Office</b>	/ Date	Office	/ Date	Office	/ Date
Case 1 of 2								
	IN	OUT	IN	OUT	IN	OUT	IN	OUT
Yag Laser Base								
Roller bar in zipper pocket								
Case 1 of 2			•					
Head/chin rest								
Case 2 of 2	IN	OUT	IN	OUT	IN	OUT	IN	OUT
Box Number 1								
Box Number 2								
Inner Black Bag #1								
Lenses	IN	OUT	IN	OUT	IN	OUT	IN	OUT
Ocular Peyman 12-12 Capsulotomy YAG								
Abraham Iridectomy Yag Laser Lens								
Drops	IN	OUT	IN	OUT	IN	OUT	IN	OUT
3 mirror gonio lens								
Eye Wash								
Tropicamide 1%								
Pilocarpine 1%								
Pilocarpine 2%								
Goniosoft								
Yag Laser Key								
Inner Black Bag #2	IN	OUT	IN	OUT	IN	OUT	IN	OUT
Power Cord								
Black Hanled 2pcs								
Track Covers 2 pcs								
YAG Laser Tower								
Optics								